



3B BlackBio Biotech India
A joint venture of Kipest India Limited, 2B BlackBio, S.L. and Bioteools B&M Labs, S.A. Madrid, Spain

**7-C, Industrial Area, Govindpura
Bhopal-462 023**

Tel. +91-755-4077847

Fax +91-755-4282659

E-mail: info@3bblackbio.com, orders@3bblackbio.com

Website: www.3bblackbio.com

3B dNTPs

Individual Deoxynucleotide Triphosphate Separate Solutions

dATP, dCTP, dGTP, dTTP

STORE AT -20°C

Description

3B dNTPs consists in separate individual solutions of the four deoxynucleotide triphosphate dATP, dACTP, dGTP, dTTP in pure water pH 7.0. The individual dNTP can be supplied at two concentrations; 10 mM or 100 mM. The dNTPs can be mixed using the four vials in a set (dATP, dCTP, dGTP, dTTP), or using only the required dNTPs (depending on the remaining dNTP to be included in the reaction as a labelled / modified form).

This individual dNTP format allows flexibility in preparation of dNTP mixtures for different applications, standard or specialised, where deoxynucleotide triphosphates are necessary for the synthesis of new DNA strands e.g. PCR, RT-PCR, DNA sequencing, DNA labelling, etc. Each dNTP can be diluted and mixed with other dNTPs to the desired user concentration.

The recommended concentration for amplification reactions is about 50-500 μM each dNTP, being the most commonly used dNTP concentration 200 μM . An increase in the concentration of dNTPs should be accompanied by an increase of the concentration of MgCl_2 , because high concentrations of dNTPs behave as potent Mg^{2+} chelating agent reducing therefore the availability of free Mg^{2+} for polymerase activity.

For amplification of DNA fragments less than 3 Kb in length do not use concentrations higher than 200 μM (each). Long amplifications, on the other hand, require higher concentrations of dNTPs, namely in the range 300-500 μM (each), as well as special DNA polymerase combinations. Each Lot of 3B dNTPs is tested under strict conditions to ensure Lot quality and Lot-to-Lot reproducibility.

For experiments using only some of the separate dNTPs vials (e.g. labeling reactions), we recommend following the procedure hereby indicated:

1. Mix in a sterile vial the required individual dNTPs.
2. Add an equimolar concentration of the labelled/modified dNTP. For diluting the dNTPs you can use sterile bidistilled water.
3. It is important that all dNTPs in a reaction mixture would be at an equimolar concentration, in order to avoid misincorporations. If the molar concentration of the labelled/modified dNTP is lower, supplement with unlabelled dNTP of the same nature (e.g. if 1 μmol of labelled dATP is added, in a solution containing 4 μmol of dCTP, 4 μmol of dGTP and 4 μmol of dTTP, add 3 μmol of unlabelled dATP).

Storage Conditions

Store at -20°C in a constant temperature freezer. Do not freeze/thaw repeatedly, for a frequent use aliquot vial content.

Order Information

dNTP Mix

PRODUCT	FORMAT	CAT N°
10 mM dNTP Mix	250 µl	3B151
10 mM dNTP Mix	2 x 250 µl	3B152
10 mM dNTP Mix	4 x 250 µl	3B153
25 mM dNTP Mix	250 µl	3B154

Individual dNTPs

PRODUCT	FORMAT	CAT N°
10 mM dATP	250 µl	3B156
10 mM dCTP	250 µl	3B157
10 mM dGTP	250 µl	3B158
10 mM dTTP	250 µl	3B159
100 mM dATP	250 µl	3B160
100 mM dCTP	250 µl	3B161
100 mM dGTP	250 µl	3B162
100 mM dTTP	250 µl	3B163

Notice to buyers/users:

Some of the applications which may be performed with this product are covered by applicable patents in certain countries. The purchase of this product does not include or provide a license to perform patented applications. Users may be required to obtain a license depending on the country and/or application.

Ed. 20 – April 11