

## WARRANTY

Products are guaranteed to conform to the quality and content indicated on each vial and external labels during their shelf life. 3B BLACKBIO BIOTECH's obligation and purchaser's rights under this warranty are limited to the replacement by 3B BLACKBIO BIOTECH of any product that is shown defective in fabrication, and that must be returned to 3B, freight prepaid, or at 3B option, replacement of the purchasing price. Any complaint on damaged goods during transport must be directed to the handling or transport agent.

Product for Research Use Only. This product must be used by qualified professionals only. It is the responsibility of the user to ascertain that a given product is adequate for a given application. Any product not fulfilling the specifications included in the product sheet will be replaced. This warranty limits our responsibility to the replacement of the product. No other warranties, of any kind, express or implied, including, without limitation, implicit warranties of commercialisation ability or adequacy for a given purpose, are provided by 3B. 3B will not be held responsible for any direct, indirect, consequential or incidental damage resulting of the use, misuses, results of the use or inability to use any product.

**Notice to users:** Some of the applications which may be performed with this product are covered by applicable patents in certain countries. The purchase of this product does not include or provide a license to perform patented applications. Users may be required to obtain a license depending on the country and/or application.

Ed.20- April 11



3B BlackBio Biotech India  
A joint venture of Kipsett India Limited, 3B BlackBio, S.L. and Biotech B&B Labs, S.A. Madrid, Spain

7-C, Industrial Area, Govindpura  
Bhopal-462 023

Tel. +91-755-4077847  
Fax +91-755-4282659

E-mail: [info@3bblackbio.com](mailto:info@3bblackbio.com),  
[orders@3bblackbio.com](mailto:orders@3bblackbio.com)

Website: [www.3bblackbio.com](http://www.3bblackbio.com)



3B BlackBio Biotech India  
A joint venture of Kipsett India Limited, 3B BlackBio, S.L. and Biotech B&B Labs, S.A. Madrid, Spain

# 3B ULTRATOOLS DNA POLYMERASE 1 U/ $\mu$ l

PRODUCT	FORMAT	REF.
UltraTools DNA Polymerase 1U/ $\mu$ l-Standard reaction Buffer	100U	3B064
UltraTools DNA Polymerase 1U/ $\mu$ l-Standard reaction Buffer	250U	3B065
UltraTools DNA Polymerase 1U/ $\mu$ l-Reaction Buffer MgCl <sub>2</sub> Free	100U	3B068
UltraTools DNA Polymerase 1U/ $\mu$ l- Reaction Buffer MgCl <sub>2</sub> Free	250U	3B069

## 1. GENERAL CONSIDERATIONS

Highly thermostable DNA polymerase. It is a recombinant, modified form of the enzyme from the thermophilic bacterium *Thermus thermophilus* expressed in *E. coli*.

3B BlackBio Biotech ULTRATOOLS DNA polymerase is suitable for applications which require a highly thermostable and processive enzyme capable of synthesising DNA strands at elevated temperatures in DNA amplification reactions or similar (e.g. primer extension), thus resolving the most complex secondary structures.

Our ULTRATOOLS DNA polymerase is also recommended for non-stringent applications (e.g. **RAPDs**). It is the enzyme of choice for applications involving **bacterial DNA sequences homologous to those found in *E. coli***.

The enzyme is free of unspecific endo- or exonucleases activities, as well as nicking activities. It does not either exhibit significant reverse-transcriptase activity. Terminal transferase activity inherent to the enzyme renders A-tailed amplification products suitable to be further used in T/A cloning approaches

*The enzyme is supplied at a concentration of 1 U/ $\mu$ l in a storage buffer. This concentration allows accurate pipetting of small amounts of the DNA polymerase, so that it is not necessary to perform ulterior dilutions.*

### Unit Definition

One unit is defined as the amount of enzyme which incorporates 10 nanomoles of dNTPs into acid-insoluble DNA within 30 minutes at 72 °C.

### Reaction buffer

Recommended reaction buffer is: 75 mM Tris HCl (pH 9.0), 2 mM MgCl<sub>2</sub> (see Note 1), 50 mM KCl, 20 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>. This reaction buffer (the so-called Standard Buffer, Ordering Information at the end) is supplied at 10X concentration together with the enzyme (either accompanying it, or included in it, in the case of aliquoted vials).

Reaction buffer can be supplied MgCl<sub>2</sub> free (the so-called Free Buffer, see Ordering Information): Mg<sup>2+</sup> ion, being the enzyme cofactor, plays a key role on polymerase activity, this is why its concentration must be optimised in certain amplification-based experiments. In this case, the MgCl<sub>2</sub> is supplied as a separate vial at 50 mM concentration. This solution must be completely thawed, vigorously vortexed and spun down in a bench-top centrifuge before use.

### Storage Conditions

Store at -20°C in a **constant temperature freezer** (i.e. do not use frost-free freezers). Under these conditions the activity of the enzyme remains unaltered over 18 months of storage. The glycerol in the storage buffer prevents freezing at -20°C.

### Reaction Conditions

After thawing the reaction buffer (and MgCl<sub>2</sub> solution, in case the "free buffer" choice is adopted), shake all vials (buffer, enzyme, 50 mM MgCl<sub>2</sub> solution) by gentle vortexing, later spin them down in a bench-top centrifuge, and eventually pipette desired volumes.

Keep all reagents on ice while they remain out of the -20°C storage freezer, otherwise enzyme activity will decrease over the time. Wear disposable gloves and make use of sterile, DNase- and RNase-free pipette tips and tubes in order to avoid contaminations and false negative results.

### Recommended enzyme volumes to be added to the reaction mix

Final reaction volumes	Recommended enzyme volumes
100 $\mu$ l	Up to 2.5 $\mu$ l
50 $\mu$ l	1-1.25 $\mu$ l
25 $\mu$ l	0.5-0.75 $\mu$ l

It is recommended to increase the enzyme units (up to three times) in order to perform certain applications such as PRINS (Primed In Situ Synthesis) or when working on long DNA fragment amplifications (longer than 2 Kb from genomic DNA).

The dNTP final concentration recommended is 200  $\mu$ M, but this figure may be altered (e.g. when unspecific amplimers occur), enlarged (e.g. long amplifications) or even unbalanced in favour of any dNTP in particular (e.g. in vitro mutagenesis experiments) depending on the intended approach. This enzyme also accepts modified dNTPs (e.g. radioactively or fluorescein labelled) as substrate.

### Notes

Note 1: at difference with the vast majority of the thermostable DNA polymerases existing in the market, our ULTRATOOLS DNA polymerase shows optimal specificity at 2 mM MgCl<sub>2</sub> final concentration (rather than 1.5 mM) in reaction buffer.